
Revision History			
Rev.	Description	Ву	Date
	Initial Release	N/A	N/A
1	Change tolerance from +/- 0.005" to 0.007" Added MicroVu inspection requirement.	N/A	09/10/08
2	cell dim060" is actually .057" & line weight 0.010" is actually 0.007", chg to meet manuf. practice	PM	08/25/11



Note: High Accuracy = \pm - 0.007" registration relative to top/left sides. Must use MicroVu for First Article and In Process Inspections.

Notes:

Cell Size: 0.060" x 0.057"
 Line Weight: 0.007"
 Font: Helvetica Bold 3 pt.

NOTICE OF PROPRIETARY RIGHTS IN DATA

THE DATA CONTAINED HEREIN IS PROPRIETARY TO GEL-PAK. THIS DATA SHALL NOT BE DUPLICATED, TRANSFERED, MADE AVAILABLE TO, OR USED BY ANY THIRD PARTY FOR ANY PURPOSE EXCEPT AS SPECIFICALLY AUTHORIZED IN WRITING BY GEL-PAK.



31398 Huntwood Ave. Hayward, CA 94544 (510) 576-2220

www.gelpak.com

Description:

GEL-PAK Print# P-0041-HA (20x20, High Accuracy)

 Part #
 P-0041-HA
 REV: 2

 Dwg #
 P-0041-HA_VR76-195.ai
 Date: 08/25/11

Scale: 1:1 Dimensions: Inches Sheet: 1 of 1